
Mölnlycke®

# Mepilex® Border Flex

WITH SAFETAC® TECHNOLOGY

Oval

STERILE EO

Single use  
Usage unique

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## Mepilex® Border Flex

Self-adherent soft silicone foam dressing

### Product description

Mepilex Border Flex is a self-adherent, absorbent dressing that maintains a moist wound environment. The waterproof outer layer protects the wound from dirt and bacteria. The dressing has a Safetac® wound contact layer that is a unique adhesive technology. It minimises pain to patients and trauma to wounds and the surrounding skin at dressing removal.

#### Mepilex Border Flex consists of:

- a wound contact layer consisting of soft silicone adhesive (Safetac) and a film carrier
- a flexible absorbent pad in three layers: a foam, a non-woven spreading layer and a layer with super absorbent fibres; the wound pad is partly perforated with Flex cut technology
- an outer film which is breathable but waterproof, providing a barrier to external contaminants

#### Dressing material content:

Silicone, polyurethane, polyacrylate, cotton, viscose, polyester and polyolefin.

### Indications for use

Mepilex Border Flex is designed for a wide range of exuding wounds such as pressure ulcers, leg and foot ulcers, traumatic wounds (e.g. skin tears) and surgical wounds. Mepilex Border Flex can also be used on dry/necrotic wounds in combination with gels.

Mepilex Border Flex reduces postoperative blistering, and may also be used as part of a prophylactic therapy to help prevent skin damage, e.g. pressure ulcers.

### Precautions

- Do not use on patients with known hypersensitivity to the ingoing materials/components of the product.
- Do not use together with oxidising agents such as hypochlorite solutions or hydrogen peroxide.
- If you see signs of infection e.g. fever or the wound or surrounding skin becoming red, warm or swollen, consult a health care professional for appropriate treatment.
- The use of dressings as part of a prophylactic therapy does not preclude the need to continue to develop and follow a comprehensive pressure ulcer prevention protocol, i.e. support surfaces, positioning, nutrition, hydration, skin care and mobility.
- Do not reuse. If reused performance of the product may deteriorate, cross contamination may occur.
- Sterile. Do not use if sterile barrier is damaged or opened prior to use. Do not re-sterilise.

### Instructions for use

Mepilex Border Flex can be used by lay persons under supervision of health care professionals.

1. Select an appropriate dressing size/shape. When applicable, choose a dressing shape optimized for the specific anatomical site.

For use on wounds:

The wound pad should overlap the dry surrounding skin by at least 1-2cm in order to protect the surrounding skin from maceration and to affix the dressing securely.

If used for prevention:

Assure that the wound pad covers the area at risk of pressure injury/tissue damage.

2. For use on wounds:

Cleanse the wound in accordance to clinical practice. Dry the surrounding skin thoroughly.

3. Remove the first release film and apply the adherent side to the wound.

4. Remove the remaining release liners and smooth down the border on the skin. Do not stretch.

Mepilex Border Flex may be left in place up to 7 days, depending on the condition of the wound and surrounding skin, or as indicated by clinical practice.

#### For use on wounds:

- Mepilex Border Flex can be used under compression bandaging.
- As Mepilex Border Flex maintains a moist wound environment, which supports debridement, there might be an initial increase in the wound size. This is normal and to be expected.
- A change in dressing regimen can result in an initial increased level of exudates, which temporarily may require an increased change frequency.

#### If used for prevention:

The area at risk of pressure ulcer/tissue damage should be inspected at regular intervals according to clinical practice.

Disposal should be handled according to local environmental procedures.

### Other information

The polyurethane foam used in the product may change colour to more yellow when it is exposed to light, air and/or heat. The colour change has no influence on product properties when used before expiry date.

If any serious incident has occurred in relation to the use of Mepilex Border Flex it should be reported to Mölnlycke Health Care.

Mepilex® and Safetac® are registered trademarks of Mölnlycke Health Care AB.

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**CANADA** Mölnlycke Health Care Inc.,  
2010 Winston Park Drive, Suite 100,  
Oakville, Ontario, L6H 5R7  
Tel: +1 905 829 1502

**CHINA** Mölnlycke Health Care  
(Shanghai) Co., Ltd.  
RM 629, No.8 Hua Jing Road,  
Wai Gao Qiao FTZ Shanghai, China  
Tel: +86 10 5128 8571

**CZECH REPUBLIC** Mölnlycke Health Care s.r.o.  
Hajkova 2747/22, 130 00 Prague 3  
Tel: +420 221 890 517 (**Reception**)  
Tel: +420 221 890 511 (**Customer Service**)

**DENMARK** Mölnlycke Health Care ApS,  
Gydevang 33, DK-3450 Allerød  
Tel: +45 48 16 82 68  
**Customer Service** Tel: +80 886 810

**ESTONIA** Mölnlycke Health Care OÜ,  
Vabaõhumuuseumitee 4-67, EE-13522 Tallin  
Tel: +372 671 1520

**FINLAND** Mölnlycke Health Care Oy,  
Pitäjänmäentie 14, 3 krs, FI-00380 Helsinki  
Tel: +358 201 622 300  
**Customer Service** Tel: 0800 113 442

**FRANCE** Mölnlycke Health Care SAS,  
13 Allée du Château Blanc  
59290 WASQUEHAL  
Tel: +33 320 122 555  
**Customer Service** Tel: 0800 910 292

**GERMANY** Mölnlycke Health Care GmbH  
Grafenberger Allee 297  
DE-40237 Düsseldorf  
Tel: +49 211 92 08 80  
**Customer Service** Tel: 0800 186 21 80

**HUNGARY** Mölnlycke Health Care Kft  
HU-1134 Budapest, Dévai u. 26-28. III. em.  
Tel: +36 1 477 3080

**INDIA** Mölnlycke Health Care India Pvt Ltd  
Universal Business Centre  
Second Floor, Universal Trade Tower  
Sector 49, Gurgaon Sohna Road  
Gurgaon 122 018, Haryana  
Tel: +91 124 4696 008

**ITALY** Mölnlycke Health Care s.r.l.,  
Via Marsala 40/C, IT-21013 Gallarate (VA)  
Tel: +39 0331 714 011

**JAPAN** Mölnlycke Health Care KK,  
4F Concieria Nishi-Shinjuku Tower's West  
6-20-7 Nishi-Shinjuku,  
Shinjuku-Ku, Tokyo 160-0023  
Tel: +81 3 6914 5004

**KOREA** Mölnlycke Health Care Korea Co., Ltd  
2nd floor, Hansol Richville.  
46 Wiryeseong-daero,  
Songpa-gu, Seoul, 05627  
South Korea  
Tel: +82 2 3789 1402

**LATVIA** Mölnlycke Health Care AB,  
Dzelzavas 120 M, Riga LV-1021  
Tel: +372 6711 522

**LITHUANIA** Mölnlycke Health Care AB,  
J. Galvydzio g. 3, LT 08236 Vilnius  
Tel: +370 5 2274 58 37

**NETHERLANDS** Mölnlycke Health Care B.V.  
Postbus 3196, 4800 DD Breda  
Tel: +31 76 521 9663

**NORWAY** Mölnlycke Health Care AS  
Postboks 6229 Etterstad, NO-0603 Oslo  
Tel: +47 22 70 63 70  
**Customer Service** Tel: 800 161 78

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Mölnlycke Health Care AG
Mölnlycke Health Care Ltd.
**USA** Mölnlycke Health Care US, LLC